

PATIENT AFTER-CARE INSTRUCTIONS



Your wound has been repaired using exofin® Precision Pen. The list below is a guide for you to understand and care for your wound following your procedure.

1 KEEP THE WOUND DRY.

You may occasionally and briefly wet your wound in the shower at the direction of your physician, but do not soak or scrub the wound area. After showering, gently blot the wound dry with a soft towel.

2 AVOID TOPICAL MEDICATIONS.

Do not apply liquid or ointment medications, lotions, creams, petroleum jelly, mineral oils or any other product to your wound while the exofin® adhesive film is in place.

3 DO NOT RUB, SCRATCH, OR PICK AT THE WOUND.

Doing so may compromise the integrity of the wound closure and cause scarring.

4 PROTECT THE WOUND FROM PROLONGED SUNLIGHT EXPOSURE.

Do not use tanning lamps while the film is in place.

5 CHECK WOUND APPEARANCE.

Some swelling, redness, and pain are common with all wounds and normally will go away as the wound heals. If swelling, redness, or pain increases, or if the wound feels warm to the touch, contact your doctor. Also contact your doctor if the wound edges reopen or separate.

6 THE ADHESIVE WILL NATURALLY SLOUGH OFF BETWEEN 7 TO 10 DAYS AFTER THE PROCEDURE.

By this time, your wound should be sufficiently healed. This information is not intended as a substitute for the advice of a physician. For more detailed information, talk to your doctor. Your doctor can choose the best treatment for you in your particular circumstances.

For additional questions and concerns, please consult your doctor.

44560054 REV2.0 V03262025