
Tourniquet Cuffs

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English

INTENDED USE

Tourniquet cuffs are used during surgical procedures to temporarily occlude the blood flow in patient's upper or lower extremities.

Location of use: Rooms suitable for medical purposes.

INDICATION / CONTRAINDICATION

Indications and contraindications depend on the intervention and therefore on the selected tourniquet cuff.

Possible indications:

- Reduction of certain fractures
- Arthroscopy of knee, hand, finger or elbow
- Bone grafts
- Kirschner wire removal
- Traumatic or nontraumatic amputations
- Tumour and cyst excisions
- Subcutaneous fasciotomy
- Nerve injuries
- Tendon repair
- Replacement or revision of joints (knee, wrist or finger)
- Correction of a hammer toe
- Podiatry

Further indications are not known.

Possible contraindications:

- Open fractures of extremities
- Post-traumatic lengthy hand reconstruction
- Severe crushing injuries
- Elbow surgery where there is excess swelling
- Severe hypertension
- Skin transplant
- Compromised circulation (e.g. peripheral artery disease)
- Diabetes mellitus

Further contraindications are not known.

In every case, the physician has to evaluate the indications and contraindications due to his or her expert knowledge.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product.
- This product must only be used by a physician, or by medically trained personnel working under the instruction of a physician.
- Prior to each use, the tourniquet cuffs must be inspected visually and a functional check must be performed (see chapter "Checks").
- The tourniquet cuffs must be reprocessed prior to each use (see chapter "Reprocessing (Cleaning, Disinfection, Sterilisation)").
- These tourniquet cuffs are designed and tested for use with tourniquet devices from the manufacturer. If the user uses tourniquet devices from other manufacturers, then the manufacturer assumes no liability for the tourniquet cuffs.
- Ensure safe connection with the device and verify continuously.

CHECKS

Visual check

Ensure that there are no cracks, breakage, etc. available on the tourniquet cuff.

Check that the hook-fastener is free from foreign body and that it adheres sufficiently.

Functional check

Perform the functional check by using the tourniquet device. Read and follow the instruction for use of the tourniquet device.

Damaged or leaky products must not be used.

USE

CAUTION

Thin padding under the cuffs is recommended to avoid pressure sores and injury to the skin.
To avoid chemical burns to the skin, seal the cuff at the distal edge. Use protective covers for leak-proof sealing and to protect the cuff from contamination.

Select proper cuff size using the colour-coded tape measure.

CAUTION

The bladder must overlap the extremity completely to ensure complete occlusion of blood flow. However, too much overlap may cause cuff rolling.

Apply the cuff tightly to the extremity (avoid wrinkling) and connect to the device.

For I.V.R.A. the minimum occlusion time must be 20 minutes.

To prevent the connection tubings from being in the operating field they should be directed proximally.

Inflate the cuff with the lowest necessary pressure (arm: max. 300 mmHg, leg: max. 400 mmHg).

CAUTION

Regarding the duration of the bloodless field, established scientific guidelines must be followed.

REPROCESSING (CLEANING, DISINFECTION, STERILISATION)

TABLE OF POSSIBLE REPROCESSING PROCEDURES

Procedures to be used	Automated procedure	Manual procedure	Sterilisation (after cleaning and disinfection)	Wipe disinfection
20-5x-xxx (max. 60 °C / 140 °F)	-	++	-	++
20-6x-xxx	++	+	++	+
20-7x-xxx (max. 60 °C / 140 °F)	-	++	-	++

++ optimal reprocessing procedure

+ permitted but suboptimal reprocessing procedure

- not permitted reprocessing procedure

GENERAL INFORMATION

- The devices are supplied non-sterile and must be reprocessed prior to first use and before each use according to the reprocessing table.
- The manufacturer is not responsible for damage resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the parameters used during every reprocessing.
- The effectiveness has been verified by an independent and accredited test laboratory.
- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the device and must be removed immediately after use.

CLEANING / DISINFECTION

Automated Cleaning / Disinfection

Notes regarding the use of washer disinfectors

- Use deionised water.

Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agent must be subsequently neutralised.
- Do not use drying agents.

Procedure validated by manufacturer:

1. The tubing connections of the tourniquet cuffs must be sealed with caps to prevent penetration of fluids into the bladder. Suitable caps are available from the manufacturer (REF 22-60-000).
2. Place devices, loosely wound, in the injector carriage.

3. Start programme with the following parameters:
 - a. Pre-rinse with deionised water at 20 °C (68 °F) and an exposure time of 1 min.
 - b. Clean at 55 °C (131 °F) with an exposure time of 5 min. using deionised water and Sekumatic® ProClean cleaning agent (dose: 0.5 % (5 ml/l)).
 - c. Neutralise with Sekumatic® FNZ at 20 °C (68 °F) with an exposure time of 2 min. (dose: 0.1 % (1 ml/l)).
 - d. Rinse with deionised water at 20 °C (68 °F) and an exposure time of 2 min.
 - e. Thermal disinfection at 93 °C (199 °F) for 5 min. with deionised water.
 - f. Dry at 100 °C (212 °F).
4. Check for visible contamination. If necessary repeat the reprocessing.
5. Remove the stoppers from the tubing connections.
6. Check the devices according to the chapter "Checks" and prepare them for sterilisation (see chapter "Sterilisation").

Manual disinfection

Notes regarding the use of disinfectants

- Prepare the disinfection solution prior to each manual cycle.
- Use deionised water.

Procedure validated by manufacturer:

1. Prepare a 2 % (30 ml/l) disinfection solution of Sekusept® Aktiv with deionised water at 20 °C (68 °F). After 15 min the disinfection solution can be used.
2. The tubing connections of the tourniquet cuffs must be sealed with suitable caps.
3. Clean the devices with soft sponges in the disinfection solution. Hard brushes and other materials which may damage the surface of the devices must not be used.
4. Place the devices into the disinfection solution and leave for 15 min. Ensure all devices are covered with solution.
5. Remove disinfection solution by rinsing thoroughly with deionised water. Remaining residues can shorten the device life span or can lead to material damage.
6. Hang the unrolled devices onto a drainer for drying. Avoid any accumulation of water.
7. Check for visible contamination. If necessary repeat the reprocessing.
8. Remove the stoppers from the tubing connections.
9. Check the devices according to the chapter "Checks" and prepare those cuffs that are suitable for sterilisation (see reprocessing table) for the sterilisation procedure (see chapter "Sterilisation").

Wipe disinfection

Wipe disinfection using conventional surface disinfectants based on alcohol (e.g., Incides N® from Ecolab Deutschland GmbH) is suitable for disinfection. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal, tuberculocidal and virucidal. After the wipe disinfection, the devices must be inspected for visible contamination. If necessary, repeat the wipe disinfection. After the wipe disinfection, the devices must be checked in accordance with the chapter "Checks".

A list of compatible cleaning agents is available at www.vbm-medical.de in the download section.

STERILISATION

Notes regarding sterilisation:

- Sterilise the tourniquet cuffs rolled loosely and fixed lightly using the straps.
- Devices or sterilisation packages must be protected from mechanical damage.

After the automated/manual cleaning and disinfection procedure, the devices must be packed into a package that is suitable for steam sterilisation.

Sterilisation validated by the manufacturer:

- The devices must be steam sterilised using a fractional vacuum process.
- Exposure time is 5 min. with a sterilisation temperature of 134 °C (273 °F).

REUSABILITY / SHELF LIFE

The end of the shelf life for reusable devices is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable devices can be reprocessed up to 50 times if used properly and if the cleaning, disinfection and sterilisation instructions are complied with. Every further reuse is the responsibility of the user (see section "Checks").

Warning

By using the devices on patients with suspected prion disease there may be a high risk of transmission. In such case the doctor must decide whether the device is either disposed of (see section "Disposal") or reprocessed according to national directives.

STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the devices must have run through the complete reprocessing process to exclude any risk for the staff of the manufacturer. For safety reasons the manufacturer preserves the right to refuse soiled or contaminated devices.

DISPOSAL

Used or damaged devices must be disposed of in accordance with the relevant national and international legal regulations.

ITEM NUMBERS

	Tourniquet Cuffs, reusable
20-54-700	Single Cuff for infant, length 20 cm / 8 inch
20-54-710	Single Cuff for child, length 30 cm / 12 inch
20-54-711	Single Cuff for arm, length 35 cm / 14 inch
20-54-712	Single Cuff for long arm, length 46 cm / 18 inch
20-54-729	Single Cuff for super long leg, length 107 cm / 42 inch
20-54-512	Single Cuff for lower leg/arm, contour shape, length 46 cm / 18 inch
20-54-522	Single Cuff for leg, contour shape, length 61 cm / 24 inch
20-54-527	Single Cuff for long leg, contour shape, length 76 cm / 30 inch
20-54-528	Single Cuff for extra long leg, contour shape, length 86 cm / 34 inch
20-64-700	Silicone Single Cuff for infant, length 20 cm / 8 inch
20-64-710	Silicone Single Cuff for child, length 30 cm / 12 inch
20-64-611	Silicone Single Cuff for arm, length 35 cm / 14 inch
20-64-612	Silicone Single Cuff for long arm, length 46 cm / 18 inch
20-64-512	Silicone Single Cuff for lower leg/arm, contour shape, length 46 cm / 18 inch
20-64-522	Silicone Single Cuff for leg, contour shape, length 61 cm / 24 inch
20-64-527	Silicone Single Cuff for long leg, contour shape, length 76 cm / 30 inch
20-64-528	Silicone Single Cuff for extra long leg, contour shape, length 86 cm / 34 inch
20-74-700	Single Slim-Cuff for infant, length 20 cm / 8 inch
20-74-710	Single Slim-Cuff for child, length 30 cm / 12 inch
20-74-711	Single Slim-Cuff for arm, length 35 cm / 14 inch
20-74-712	Single Slim-Cuff for long arm, length 46 cm / 18 inch
20-74-715	Single Slim-Cuff for lower leg/arm, contour shape, length 46 cm / 18 inch
20-74-722	Single Slim-Cuff for leg, contour shape, length 61 cm / 24 inch
20-74-727	Single Slim-Cuff for long leg, contour shape, length 76 cm / 30 inch
20-74-728	Single Slim-Cuff for extra long leg, contour shape, length 86 cm / 34 inch
20-74-729	Single Slim-Cuff for super long leg, contour shape, length 107 cm / 42 inch

20-50-700	Double Cuff for infant, length 20 cm / 8 inch
20-50-710	Double Cuff for child, length 30 cm / 12 inch
20-50-711	Double Cuff for arm, length 5 cm / 14 inch
20-50-712	Double Cuff for long arm, length 46 cm / 18 inch
20-50-722	Double Cuff for leg, length 61 cm / 24 inch
20-50-727	Double Cuff for long leg, length 76 cm / 30 inch
20-50-728	Double Cuff for extra long leg, length 86 cm / 34 inch
20-50-729	Double Cuff for super long leg, length 107 cm / 42 inch
20-60-711	Silicone Double Cuff for arm, length 35 cm / 14 inch
20-60-712	Silicone Double Cuff for long arm, length 46 cm / 18 inch
20-60-722	Silicone Double Cuff for leg, length 61 cm / 24 inch
	Accessories
22-60-000	Sealing Cap for Tourniquet Cuffs
20-11-111	Cuff Cover for Tourniquet Cuff small, non sterile
20-11-111-1	Cuff Cover for Tourniquet Cuff small, sterile
20-11-333	Cuff Cover for Tourniquet Cuff medium, non sterile
20-11-333-1	Cuff Cover for Tourniquet Cuff medium, sterile
20-11-222	Cuff Cover for Tourniquet Cuff large, non sterile
20-11-222-1	Cuff Cover for Tourniquet Cuff large, sterile

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