




Roll-On Cuff

DE	ROLLMANSCHETTE Gebrauchsanweisung 4 - 6	NL	ROLMANCHET Gebruiksaanwijzing 22 - 24
EN	ROLL-ON CUFF Instructions for use 7 - 9	RU	ЭЛАСТИЧНАЯ МАНЖЕТА Инструкция по применению ... 25 - 28
DA	RULLEMANCHET Bruganvisning 10 - 12	SV	RULLMANSCHETT Bruksanvisning 29 - 31
ES	MANGUITO "ROLL-ON" Instrucciones de uso 13 - 15	TR	RULO KAF Kullanım kılavuzu 32 - 34
FR	BRASSARD ROULANT Manuel d'utilisation 16 - 18		SYMBOL DESCRIPTION 35 - 36
IT	BRACCIALE ARROTOLABILE Istruzioni per l'uso 19 - 21		

Deutsch

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English

INTENDED USE

The Roll-On Cuff is used to exsanguinate a patient's upper or lower extremity before using a Tourniquet Cuff.

Clinical benefit: the Roll-On Cuff creates a bloodless field in combination with a Tourniquet Cuff.

Patient target group: patients who require surgery at the upper or lower extremities.

Place of use: rooms which are suitable for surgical interventions.

INDICATIONS

- Surgeries in bloodless field
- No other indications are known.

CONTRAINDICATIONS

- Open fractures of the extremities
- Severe crushing
- Severe hypertension
- Skin transplant - impaired circulation (e.g. peripheral arterial disease)

No other contraindications are known.

SAFETY INSTRUCTIONS

- Read and follow the instructions for use carefully before using the product and keep them for future reference.
- This product must only be used by medically trained personnel.
- The user and/or patient must report all serious adverse events that occurred in connection with the product to the manufacturer and competent authorities of the EU member state (or report to the competent authorities of the country if an event occurs outside of the EU) in which the user and/or patient is located.
- Before every use, the product must be inspected visually for damage (cracks, breakage, etc.) (see "Checks"). A faulty product must not be used.
- It is not permitted to make any changes to the product.
- The Roll-On Cuff is designed and tested for use with the hand inflator (not included in the scope of delivery) from VBM Medizintechnik GmbH. If the user uses a hand inflator from another manufacturer, VBM Medizintechnik GmbH shall not assume any liability for the product.
- The reusable product is supplied non-sterile and must be cleaned, disinfected and sterilised prior to every use.
- The product is not suitable for MRI.

PREPARATION PRIOR TO USE

- ▶ Choose the correct size of Roll-On Cuff (a) using the colour-coded measuring tape (b) for the respective extremity (see "Product specifications", Figure 1).

CAUTION

The maximum given circumference of the extremity must not be exceeded.

APPLICATION

CAUTION

- The Roll-On Cuff must not be used as a tourniquet.
- It is essential to insert the needle of the connecting tubing (e) vertically into the valve of the Roll-On Cuff. Insertion at an angle may damage the valve irreparably and thus cause it to leak.
- Never change the pressure during application.
- Ensure that the Roll-On Cuff is never inflated to more than 130 mmHg, otherwise it will be overstretched. The material could be irreparably damaged.

- ▶ Lubricate the needle (d) with water before inserting.
- ▶ Insert the needle (d) vertically into the valve (c) (black dot) (Figure 2).
- ▶ Inflate the Roll-On Cuff (a) with the hand inflator (h) to 130 mmHg (green mark on the scale (f)).
- ▶ Separate the needle (d) from the valve (c).
- ▶ Roll the Roll-On Cuff cuff over the extremity up to the deflated Tourniquet Cuff to achieve an exsanguinated field (Figure 3).

Removing the Roll-On Cuff

CAUTION

Never deflate the Roll-On Cuff before unrolling, otherwise it cannot be unrolled.

- ▶ After inflating the Tourniquet Cuff, unroll the Roll-On Cuff from the extremity (Figure 4).
- ▶ After use, release the pressure in the Roll-On Cuff with the release valve (g) of the hand inflator (h).

REPROCESSING - ROLL-ON CUFF (CLEANING, DISINFECTION, STERILISATION)

GENERAL INFORMATION

- The manufacturer is not responsible for damage resulting from improper reprocessing.
- The user is obliged to validate their procedures or the devices and accessories and to comply with the validated parameters during every reprocessing.
- It is recommended to use an automated procedure due to a significantly higher efficiency of cleaning and disinfection.
- The effectiveness has been verified by an independent and accredited test laboratory.

- To achieve effective reprocessing, heavy soiling must not be allowed to dry on the product and must be removed immediately after use.
- Silicone products must not come into contact with oil and grease.
- The measuring tape must be disinfected according to "Wipe disinfection - measuring tape" (must not be sterilised!).

CLEANING / DISINFECTION

Automated cleaning / disinfection

Notes regarding the use of washer disinfectors

- ▶ Use de-ionised water.
 - ▶ Use a washer disinfectant which complies with ISO 15883.
- A Miele G7882 with Cabinet E435/3 was used for validation.

Notes regarding the use of cleaning agents and disinfectants

- When alkaline cleaning agents are used, the agent must be subsequently neutralised.
- Do not use drying agents.

Procedure validated by the manufacturer:

1. The air must be pressed out of the Roll-On Cuff via the release valve on the hand inflator.
2. Position the products in the injector carriage:

Product / Accessories	Positioning in the injector carriage
• Roll-On Cuff • Connecting tubing of the hand inflator	Injector jet

CAUTION

Ensure effective rinsing of all hollow spaces.

3. Start program with the following parameters:
 - a. Pre-rinse with de-ionised water at 20 °C (68 °F) and an exposure time of 1 min.
 - b. Cleaning at 55 °C (131 °F) with an exposure time of 5 min using de-ionised water and "Sekumatic® Pro-Clean" cleaning agent (dose: 0.5 % (5 ml/l)).
 - c. Neutralisation with Sekumatic® FNZ at 20 °C (68 °F) with an exposure time of 2 min (dose: 0.1 % (1 ml/l)).
 - d. Rinse with de-ionised water (max. 100 CFU/ml) at 20 °C (68 °F) and an exposure time of 2 min.
 - e. Thermal disinfection at 93 °C (199 °F) for 5 min with de-ionised water.
 - f. Dry at 100 °C (212 °F).
4. Check for visible contamination. If necessary, repeat reprocessing.
5. Check the products according to the section "Checks" and prepare them for sterilisation (see "Packaging").

Manual cleaning / disinfection

Notes regarding the use of cleaning agents and disinfectants

- Prepare the cleaning and disinfection solutions prior to each manual cycle.
- Use de-ionised water.

Procedure validated by the manufacturer:

1. Prepare a 2 % (30 ml/l) cleaning and disinfection solution of Sekusept® Aktiv with de-ionised water at 20 °C (68 °F). Stir several times for the next 15 min. The cleaning and disinfection solution is ready for use after 15 min.
2. The air must be pressed out of the Roll-On Cuff via the release valve on the hand inflator.
3. Clean the product with soft sponges in the cleaning and disinfection solution, the product must be fully immersed for at least 20 s. Hard brushes and other materials which may damage the silicone surface must not be used.
4. Place the product into the cleaning and disinfection solution and leave for 15 min.
5. Remove cleaning and disinfection solution by rinsing thoroughly with de-ionised water (max. 100 CFU/ml). Remaining residues can shorten the product shelf life or can lead to material damage.
6. Dry the product. Avoid any accumulation of water.
7. Check for visible contamination. If necessary, repeat reprocessing.
8. Check the products according to the section "Checks" and prepare them for sterilisation (see "Packaging").

CHECKS

After disinfection, the product must be inspected for damage (cracks, breakage, etc.). A faulty product must be disposed of.

PACKAGING

- ▶ The air must be pressed out of the Roll-On Cuff via the release valve on the hand inflator.
- ▶ After the cleaning and disinfection procedure, the product must be packed into a sterile barrier system suitable for steam sterilisation. The sterile barrier systems must comply with ISO 11607-1. A sterile barrier system with single packaging was used for validation.

STERILISATION

Notes regarding sterilisation:

- Product or sterilisation packaging must be protected from mechanical damage.
- Use a steriliser which complies with requirements of DIN EN 285 or DIN EN 13060.

Sterilisation validated by the manufacturer

- The product must be steam sterilised by a fractionated vacuum process.
- Exposure time is 5 min with a sterilisation temperature of 134 °C (273 °F).
- The drying time is 10 min.
- A Systec HX-320 autoclave was used for validation.

REUSABILITY

The end of the shelf life for reusable products is basically determined by wear and damage caused by its use. Considering the maximum shelf life of 5 years, reusable products can be reprocessed up to 50 times if used properly and if the cleaning, disinfection and sterilisation instructions are complied with.

Every further reuse is the responsibility of the user (see section "Checks").

WARNING

By using the products on patients with suspected prion disease there may be a high risk of transmission. In such case the physician decides if the product is either disposed of (see section "Disposal") or reprocessed according to national directives.

WIPE DISINFECTION - MEASURING TAPE

CAUTION

- The product must not be mechanically or manually reprocessed or sterilised.
- Do not immerse the product in liquids.
- In case of unsuccessful wipe disinfection or heavy contamination, the product must be disposed of as specified in "Disposal".

Wipe disinfection must be performed with commercially available surface disinfectants based on alcohol. When choosing a product for disinfection, a disinfectant with an appropriate spectrum of action must be used: bactericidal, yeasticidal and virucidal. After wipe disinfection, inspect the product for visible contamination. If necessary, repeat wipe disinfection. After the wipe disinfection, check the product (see "Checks").

SHELF LIFE

The shelf life of the product is 5 years and it can be reprocessed 50 times within the shelf life.

STORAGE AND TRANSPORT CONDITIONS

CAUTION

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

SERVICE

Prior to the return of medical devices for complaints/repairs, the products must have run through the complete reprocessing process to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer reserves the right to refuse soiled or contaminated products.

DISPOSAL

The used or faulty product must be disposed of in accordance with the applicable national and international legal regulations.

PRODUCT SPECIFICATIONS

Packaging unit: 1

REF	#	Colour coding	Circumference of extremity	Weight (incl. packaging)
21-95-711	1	white	20 - 28 cm (8 - 11 inch)	144 g (5 oz)
21-95-712	2	light blue	26 - 33 cm (10 - 13 inch)	277 g (9.8 oz)
21-95-719	3	yellow	31 - 40 cm (12 - 16 inch)	295 g (10.4 oz)
21-95-722	4	red	38 - 49 cm (15 - 19 inch)	338 g (11.9 oz)
21-95-727	5	green	47 - 60 cm (18 - 24 inch)	373 g (13.2 oz)
21-95-729	6	dark blue	58 - 70 cm (23 - 28 inch)	595 g (21.0 oz)

SCOPE OF DELIVERY

- Roll-On Cuff
- Measuring tape

MATERIAL DATA

Roll-On Cuff	Cuff: silicone Support ring: aluminium
Measuring tape	Synthetic fibre paper

ACCESSORIES

REF 21-95-300	Hand inflator with pressure gauge
---------------	-----------------------------------

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