

Designed to Optimize Wound Closure

- Tinted purple for visualization during application
- Contoured tip for better application control
- An effective microbial barrier



skin **affix**[™]



Strong and Lasting Wound Closure

skin affix is a topical skin adhesive intended to hold closed easily approximated skin edges of wounds from:

- Surgical incisions
- Minimally invasive surgical procedures
- Simple, thoroughly cleansed trauma-induced lacerations

Performance


- Sterile 2-octyl cyanoacrylate
- High viscosity formulation
- Effective microbial barrier protection
- Fast setting
- Naturally sloughing (5-10 days)

Application


- Versatile – foam tip allows for broad and narrow strokes
- Promoting controlled delivery – transparent applicator
- Easy to use – no assembly required

Patient Comfort

- Low exothermic reaction*
- Flexible and water-resistant



Before



After one layer application



Item No.	Description	Size	Pkg
M1204Z	skin affix topical skin adhesive	0.4 ml	12/box
M1204	skin affix topical skin adhesive	0.4 ml	72/case



SourceMark LLC