



TO PLACE AN ORDER:

866-528-7001 phone

866-528-6008 fax

customerservice@sourcemarkusa.com

Your needs...our focus

Standard Disposable Laryngeal Mask Airway



- No aperture bars for easier insertion of diagnostic tools
- Anatomically correct curve
- Sturdy tube for ease of insertion
- Available in full range of sizes, neonate through adult
- All seven sizes at the same convenient price
- Easy product transition

Silicone Disposable Laryngeal Mask Airway



- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, neonate through adult
- All seven sizes at the same convenient price
- Disposable designation by green cuff color
- No aperture bars for easier insertion of diagnostic tools
- Supple cuff material for less irritation
- QuikPak available with syringe & lubricant

- FDA listed and registered –
- CE and ISO Quality System certified –
- Clinically tested and proven –
- Most cost effective laryngeal mask in the market –
- Individual Sterile Package, 10/box –



TO PLACE AN ORDER:



866-528-7001 phone

866-528-6008 fax

customerservice@sourcemarkusa.com

Standard/ Silicone Disposable Laryngeal Mask Airway

MASK SIZE	1	1 ^{1/2}	2	2 ^{1/2}	3	4	5
CATALOG							
Standard	M0343	M0353	M0363	M0373	M0313	M0323	M0333
Silicone	M0342	M0352	M0362	M0372	M0312	M0322	M0332
QuikPak*	MX0342	MX0352	MX0362	MX0372	MX0312	MX0322	MX0332
PATIENT SELECTION GUIDELINES	Neonates/ Infants up to 5 kg	Infants 5-10 kg	Infants/ Children 10-20 kg	Children 20-30 kg	Children 30-50 kg	Adults 50-70 kg	Adults 70-100 kg
MAXIMUM CUFF INFLATION VOLUME (AIR)	up to 4 ml	up to 7 ml	up to 10 ml	up to 14 ml	up to 20 ml	up to 30 ml	up to 40 ml
MAXIMUM FOB	2.7 mm	3.0 mm	3.5 mm	4.0 mm	5.0 mm	5.0 mm	5.0 mm
MAXIMUM ETT	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm Cuffed	6.0 mm Cuffed	7.0 mm Cuffed

*QuikPak is SourceMark's MX series for fast and convenient usage. It includes one syringe  and one pack of 5.0g (.17oz) of lubricating jelly. 

Our laryngeal mask airways have been certified by the CE and ISO Quality Systems. ISO 13485