



## TO PLACE AN ORDER:

**866-528-7001** phone

**866-528-6008** fax

[customerservice@sourcemarkusa.com](mailto:customerservice@sourcemarkusa.com)

## Your needs...our focus

The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

### Standard High-Low Cuffed ET Tube



### Standard Un-cuffed ET Tube



### Oral Cuffed Preformed ET Tube



### Oral Un-cuffed Preformed ET Tube



### Nasal Cuffed Preformed ET Tube



### Nasal Un-cuffed Preformed ET Tube



## CUFFED:

- High volume low-pressure cylindrical cuff distributes even pressure on the tracheal wall reducing chance of tissue damage
- Large cuff surface provides better seal
- Durable construction to ensure the reliability of the tube and cuff during procedures

## TUBE:

- Thermosmoothing of the tube tips and Murphy eye to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

## GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized

**TO PLACE AN ORDER:**

866-528-7001 phone

866-528-6008 fax

customerservice@sourcemarkusa.com

**Standard High-Low Cuffed ET Tube** (10 ea per box / 10 boxes per case)

| TUBE SIZE | 4.0    | 4.5    | 5.0    | 5.5    | 6.0    | 6.5    | 7.0    | 7.5    | 8.0    | 8.5    | 9.0    |
|-----------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| CATALOG   | M0440C | M0445C | M0450C | M0455C | M0460C | M0465C | M0470C | M0475C | M0480C | M0485C | M0490C |

**Standard Un-cuffed ET Tube** (10 ea per box / 10 boxes per case)

| TUBE SIZE | 2.0    | 2.5    | 3.0    | 3.5    | 4.0    | 4.5    | 5.0    |
|-----------|--------|--------|--------|--------|--------|--------|--------|
| CATALOG   | M0420U | M0425U | M0430U | M0435U | M0440U | M0445U | M0450U |

**Oral Cuffed Preformed ET Tube** (10 ea per box / 10 boxes per case)

| TUBE SIZE | 3.0    | 3.5    | 4.0    | 4.5    | 5.0    | 5.5    | 6.0    | 6.5    | 7.0    | 7.5    | 8.0    | 8.5    | 9.0 |
|-----------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|-----|
| CATALOG   | M0530C | M0535C | M0540C | M0545C | M0550C | M0555C | M0560C | M0565C | M0570C | M0575C | M0580C | M0585C | NA  |

**Oral Un-cuffed Preformed ET Tube** (10 ea per box / 10 boxes per case)

| TUBE SIZE | 3.0    | 3.5    | 4.0    | 4.5    | 5.0    | 5.5    |
|-----------|--------|--------|--------|--------|--------|--------|
| CATALOG   | M0530U | M0535U | M0540U | M0545U | M0550U | M0555U |

**Nasal Cuffed Preformed ET Tube** (10 ea per box / 10 boxes per case)

| TUBE SIZE | 6.0    | 6.5    | 7.0    | 7.5    | 8.0    | 8.5 | 9.0 |
|-----------|--------|--------|--------|--------|--------|-----|-----|
| CATALOG   | M0660C | M0665C | M0670C | M0675C | M0680C | NA  | NA  |

**Nasal Un-cuffed Preformed ET Tube** (10 ea per box / 10 boxes per case)

| TUBE SIZE | 3.0    | 3.5    | 4.0    | 4.5    | 5.0    | 5.5    | 6.0    | 6.5    | 7.0    |
|-----------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| CATALOG   | M0630U | M0635U | M0640U | M0645U | M0650U | M0655U | M0660U | M0665U | M0670U |

**Aluminum Stylet**

(20 ea per box / 20 boxes per case)

| STYLET SIZE | 6FR   | 10FR  | 14FR  |
|-------------|-------|-------|-------|
| CATALOG     | M0906 | M0910 | M0914 |

**Indications for usage:**

- Airway management during procedures
- Anesthesia

