



## TO PLACE AN ORDER:

866-528-7001 phone

866-528-6008 fax

[customerservice@sourcemarkusa.com](mailto:customerservice@sourcemarkusa.com)

Your needs...our focus

## Standard Disposable Laryngeal Mask Airway



- No aperture bars for easier insertion of diagnostic tools
- Anatomically correct curve
- Sturdy tube for ease of insertion
- Available in full range of sizes, neonate through adult
- All seven sizes at the same convenient price
- Easy product transition

## Silicone Disposable Laryngeal Mask Airway



- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, neonate through adult
- All seven sizes at the same convenient price
- Disposable designation by green cuff color
- No aperture bars for easier insertion of diagnostic tools
- Supple cuff material for less irritation
- QuikPak available with syringe & lubricant

- FDA listed and registered –
- CE and ISO Quality System certified –
- Clinically tested and proven –
- Most cost effective laryngeal mask in the market –
- Individual Sterile Package, 10/box –



SourceMark<sub>LLC</sub>

**TO PLACE AN ORDER:**



**866-528-7001** phone

**866-528-6008** fax

[customerservice@sourcemarkusa.com](mailto:customerservice@sourcemarkusa.com)

# Standard/ Silicone Disposable Laryngeal Mask Airway

MASK SIZE	1	1 <sup>1/2</sup>	2	2 <sup>1/2</sup>	3	4	5
<b>CATALOG</b>							
<b>Standard</b>	<b>M0343</b>	<b>M0353</b>	<b>M0363</b>	<b>M0373</b>	<b>M0313</b>	<b>M0323</b>	<b>M0333</b>
<b>Silicone</b>	<b>M0342</b>	<b>M0352</b>	<b>M0362</b>	<b>M0372</b>	<b>M0312</b>	<b>M0322</b>	<b>M0332</b>
<b>QuikPak*</b>	<b>MX0342</b>	<b>MX0352</b>	<b>MX0362</b>	<b>MX0372</b>	<b>MX0312</b>	<b>MX0322</b>	<b>MX0332</b>
<b>PATIENT SELECTION GUIDELINES</b>	Neonates/ Infants up to 5 kg	Infants 5-10 kg	Infants/ Children 10-20 kg	Children 20-30 kg	Children 30-50 kg	Adults 50-70 kg	Adults 70-100 kg
<b>MAXIMUM CUFF INFLATION VOLUME (AIR)</b>	up to 4 ml	up to 7 ml	up to 10 ml	up to 14 ml	up to 20 ml	up to 30 ml	up to 40 ml
<b>MAXIMUM FOB</b>	2.7 mm	3.0 mm	3.5 mm	4.0 mm	5.0 mm	5.0 mm	5.0 mm
<b>MAXIMUM ETT</b>	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm Cuffed	6.0 mm Cuffed	7.0 mm Cuffed

\*QuikPak is SourceMark's MX series for fast and convenient usage. It includes one syringe  and one pack of 5.0g (.17oz) of lubricating jelly. 

Our laryngeal mask airways have been certified by the CE and ISO Quality Systems. ISO 13485