


LiquiBand® Rapid™
Topical skin adhesive




30 is the
new 90

to faster wound
closure and greater
efficiency in the OR



Discover LiquiBand® Rapid™—the new topical skin adhesive that dries in **as fast as 30 seconds**¹ while other products dry in 90 seconds or more²



LIQUIBAND RAPID

OTHER PRODUCTS



Now more than ever, working smarter and faster in the OR is critical. LiquiBand Rapid with Accelerated Drying Technology answers the call.



Spend less time waiting for the adhesive to dry.



Improve OR turnover times.



Care for more patients, more effectively.

Designed with the surgeon's needs in mind

User-centered design features help you close the wound with ease, consistency and accuracy no matter the setting.




POROUS FOAM TIP

helps control the application of the adhesive without dripping.



ACCELERANT-COATED TIP

speeds up the adhesive drying process and helps reduce patient turnover time.



2-OCTYL CYANOACRYLATE

formulation makes a stronger, flexible film for more durable, long-lasting wound closure.⁴



HIGHLIGHTER SHAPE


allows the user to apply broad or thin lines more accurately.



WINGED APPLICATOR

protects the user's fingers from the risk of glass shards and makes activating the adhesive safe and easy to control with one hand⁵


It's more than fast. LiquiBand Rapid provides the strong microbial barrier you expect with 2 Octyl cyanoacrylates. Note: in vitro studies show effective microbial barrier with no detectable microorganism penetration.³



Why LiquiBand Rapid is the stronger choice

Greater wound closure strength.⁶


Stronger and more durable to help prevent dehiscence.



Follows ASTM F2458-05 standard test method. Two substrates were placed opposing each other and bonded together with adhesive at the joint line to simulate a wound closure. Force was applied to separate the substrates and measured.

Better lap shear strength.⁶


Greater flexibility on the skin for long-lasting wound closure and maintenance of the microbial barrier.



Follows ASTM F2255-05 standard test method by bonding two overlapping surfaces with the adhesive. Force was applied to pull apart the surfaces and measured.

Easier to use, safe to activate.⁶

Requires less force to activate and protects the user's fingers throughout application.



Force applied to activate the device was mechanically measured, simulating the amount of pressure a clinician would need to apply.

Take the 30 is the new 90 challenge

See for yourself how LiquiBand Rapid dries in as fast as 30 seconds while other products dry in 90 seconds or more.



Item No.	Description	Pkg.
M1207	Liquiband Rapid topical skin adhesive, 0.8 mL	36/cs

References: 1. Test data on file. 2. "The-DERMABOND-Family-of-Topical-Skin-Adhesives" information sheet (039696_171101) 3-6. Test data on file. LiquiBand, LiquiBand Rapid and logos are trademarks and registered trademarks of Advanced Medical Solutions Limited.

© 2020 All rights reserved. Sourcemark is a registered trademark of Sourcemark, LLC. MKT19W1055260 / e20477 / 3
Distributed by SourceMark, LLC | 302 Innovation Drive, Suite 410, Franklin, TN 37067 | 866-528-7001 | www.sourcemarkusa.com

