





Tourniquet Accessories



Measuring Tape


For selection of VBM tourniquet cuffs

- Color coded, for easy identification of cuff selection

Order information

Measuring Tape / For single use

REF	Box
20-95-222	50



Tourniquet Sleeve


For patient protection

- Protects soft tissues underneath the cuff
- Reduces wrinkling, pinching and shearing of soft tissues

Order information

Tourniquet Sleeve / For single use

Extremity	Color	Clean packed	Sterile packed	Box
≤20 cm (8")	○	REF 23-00-700	REF 23-00-700-1	10
≤30 cm (12")	●	REF 23-00-710	REF 23-00-710-1	10
≤35 cm (14")	●	REF 23-00-711	REF 23-00-711-1	10
≤46 cm (18")	●	REF 23-00-712	REF 23-00-712-1	10
≤61 cm (24")	●	REF 23-00-722	REF 23-00-722-1	10
≤76 cm (30")	●	REF 23-00-727	REF 23-00-727-1	10
≤86 cm (34")	●	REF 23-00-728	REF 23-00-728-1	10
≤107 cm (42")	●	REF 23-00-729	REF 23-00-729-1	10



Tourniquet Drape


For patient protection

- Isolates the cuff from the sterile field
- Seals against fluid accumulation underneath the cuff to protect the patient from burns associated with skin preparation

Order information

Tourniquet Drape / For single use

Extremity	Clean packed	Box	Sterile packed	Box
≤45 cm (≤18")	REF 24-00-100	40	REF 24-00-100-1	10
35 - 76 cm (14" - 30")	REF 24-00-200	40	REF 24-00-200-1	10
61 - 107 cm (24" - 42")	REF 24-00-300	40	REF 24-00-300-1	10



Tourniquet Cuff Cover


For cuff protection

- Keeps the cuff free from soiling

Order information

Tourniquet Cuff Cover / For single use


Extremity	Clean packed	Sterile packed	Box
≤42 cm (≤16 1/2")	REF 20-11-111	REF 20-11-111-1	50
40 - 58 cm (16" - 23")	REF 20-11-333	REF 20-11-333-1	50
52 - 70 cm (20" - 28")	REF 20-11-222	REF 20-11-222-1	50



Esmarch Bandages

For exsanguination

- Available in 4 sizes
- Made of synthetic rubber and powder-free (talcum-free) to avoid possible allergic reactions
- Textured surface for a good grip
- Elastic properties for easy wrapping
- Robust and tear resistant



Order information

Esmarch Bandage / For single use, sterile packed, color: light blue

Width	Length	REF	Box
6 cm (2 3/8")	365 cm (12 ft)	20-21-060-1	10
8 cm (3 1/8")	365 cm (12 ft)	20-21-080-1	10
10 cm (4")	365 cm (12 ft)	20-21-100-1	10
15 cm (6")	365 cm (12 ft)	20-21-150-1	10

Esmarch Bandage / Reusable, autoclavable, color: dark blue

Width	Length	REF	Box
6 cm (2 3/8")	500 cm (16 ft)	20-20-060	5
8 cm (3 1/8")	500 cm (16 ft)	20-20-080	5
10 cm (4")	500 cm (16 ft)	20-20-100	5
15 cm (6")	500 cm (16 ft)	20-20-150	5



Roll-On Cuff


For exsanguination

- Allow fast and easy exsanguination within seconds
- Inflatable exsanguination system
- 100% silicone medical grade
- Quick and efficient procedure
- Full range of sizes available

Order information

Roll-On Cuff / Reusable

Size	Extremity	Color	REF	Box
2	26 - 33 cm (10" - 13")	Blue	21-95-712	1
3	31 - 40 cm (12" - 16")	Yellow	21-95-719	1
4	38 - 49 cm (15" - 19")	Red	21-95-722	1
5	47 - 60 cm (18" - 24")	Green	21-95-727	1
6	58 - 70 cm (23" - 28")	Dark Blue	21-95-729	1



Additional information



Application Video
Tourniquet Sleeve



Application Video
Tourniquet Drape



Application Video
Tourniquet Cuff
Cover



Application Video
Esmarch Bandages



Application Video
Roll-On Cuff

Accessories for Roll-On Cuff



Hand Inflator for Roll-On Cuff / Manometer 200 mmHg, Smooth Connecting Tubing (100 cm (40")) with positive locking connectors and inflation needle

REF	Box
21-95-300	1




Smooth Connecting Tubing / With positiv locking connectors and inflation needle, (100 cm (40"))

REF	Box
21-95-302	1



Brake-Block / Reusable, material: Silicone

Size	REF	Box
Small	21-95-200	1
Large	21-95-201	1



Measuring Tape for VBM Roll-On Cuff / Color coded, for easy identification of Roll-On Cuff

REF	Box
21-95-111	10

Order information



302 Innovation Drive | Suite 410

Franklin, TN 37067

USA

Phone: 866.528.7001

Fax: 866.528.6008

e-mail: customerservice@sourcemarkusa.com

www.sourcemarkusa.com