

Norm-O-Temp

NORMOTHERMIA SYSTEM

The Norm-O-Temp[®] is a whole body system that allows conductive warming therapy to be administered. It can be used in the operating room, pre-op, recovery or emergency department and when paired with Gelli-Roll[®] makes the perfect Temperature Management device for your OR.

- Offers an underbody solution that allows full access to the patient
- Cost savings solution for managing normothermia when combined with the Gelli-Roll
- Ability to pre-condition the water helps save time for the caregiver
- Small footprint allows for easy placement of the device



Norm-O-Temp®

NORMOTHERMIA SYSTEM

Ordering Information

Cat No.	Description	Quantity
111Z	Norm-O-Temp®	115v
111W	Norm-O-Temp®	230v
118	Stand with I.V. Pole	ea.
119	Low Profile Stand	ea.
Reusable Connecting Hoses (Must purchase 286 hoses separately to use blankets with Norm-O-Temp®)		
286	9' Reusable Connecting Hose (2.7m)	ea.
286-18	18' Reusable Connecting Hose (5.5m)	ea.
286-27	27' Reusable Connecting Hose (8.2m)	ea.
Gelli-Roll® (Must purchase 286 hose separately to use Gelli-Roll® with Norm-O-Temp®)		
195P	Gelli-Roll® - Adult 73.5in x 21.5in x .625in (186.7cm x 54.6cm x 1.6cm)	ea.
194P	Gelli-Roll® - Pediatric 31.25in x 23.75in x .625in (79.4cm x 60.3cm x 1.6cm)	ea.
193P	Gelli-Roll® - Infant 20in x 13in x .625in (50.8cm x 33cm x 1.6cm)	ea.

Product Specifications:

Dimensions:

9in x 15in x 18.5in
(22.9cm x 38.1cm x 46.9cm)

Weight:

Empty: 33.5 lbs. (15.2 kg)
Filled: 45.5 lbs. (20.6 kg)

Water Temperature Range:

32°C - 42°C
(89.6°F - 107.6°F)
(Heating Only)

Controller Accuracy:

Water Temp: $\pm 0.6^{\circ}\text{C}$ ($\pm 1.08^{\circ}\text{F}$)
Water Temperature Display

Range:

0°C - 52°C (32°F - 126°F)

Safety Alarms:

Audible & Visible alarms for
over temperature limits.



302 Innovation Drive, Suite 410
Franklin, TN 37067 | 866.528.7001
www.sourcemarkusa.com