




ANESTHESIA / AIRWAY MANAGEMENT



Silicone Disposable Laryngeal Mask Airway



FEATURES AND BENEFITS

- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, neonate through adult
- Color coded pilot balloon for quick size identification
- No aperture bars for easier insertion of diagnostic tools
- Supple cuff material for less irritation

MASK SIZE	DESCRIPTION	ITEM #	QTY
1	Neonates/Infants up to 5 kg	M0342	10/box
1½	Infants 5-10 kg	M0352	10/box
2	Infants/Children 10-20 kg	M0362	10/box
2½	Children 20-30 kg	M0372	10/box
3	Children 30-50 kg	M0312	10/box
4	Adults 50-70 kg	M0322	10/box
5	Adults 70-100 kg	M0332	10/box

PEDIATRIC

Standard Disposable Laryngeal Mask Airway



FEATURES AND BENEFITS

- No aperture bars for easier insertion of diagnostic tools
- Anatomically correct curve
- Sturdy tube for easier insertion
- Easy product transition
- Available in full range of sizes, neonate through adult

MASK SIZE	DESCRIPTION	ITEM #	QTY
1	Neonates/Infants up to 5 kg	M0343	10/box
1½	Infants 5-10 kg	M0353	10/box
2	Infants/Children 10-20 kg	M0363	10/box
2½	Children 20-30 kg	M0373	10/box
3	Children 30-50 kg	M0313	10/box
4	Adults 50-70 kg	M0323	10/box
5	Adults 70-100 kg	M0333	10/box

PEDIATRIC



Quikpak Silicone Disposable Laryngeal Mask Airway



FEATURES AND BENEFITS

- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, infants through adult
- Color coded pilot balloon for quick size identification
- No aperture bars for easier insertion of diagnostic tools
- Supple cuff material for less irritation

MASK SIZE	DESCRIPTION	ITEM #	QTY
1½	Infants 5-10 kg	MX0352	10/box
2	Infants/Children 10-20 kg	MX0362	10/box
2½	Children 20-30 kg	MX0372	10/box
3	Children 30-50 kg	MX0312	10/box
4	Adults 50-70 kg	MX0322	10/box
5	Adults 70-100 kg	MX0332	10/box

PEDIATRIC

*QuikPak is SourceMark's MX series for fast and convenient usage.
It includes one syringe and one pack of 5.0g (.17oz) of lubricating jelly.

Reusable Laryngeal Mask Airway



FEATURES AND BENEFITS


- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, neonate through adult
- Multi-use convenience
- No aperture bars for easier insertion of diagnostic tools and cleaning
- Supple cuff material for less irritation

MASK SIZE	DESCRIPTION	ITEM #	QTY
1	Neonates/Infants up to 5 kg	M0341	1 each
1½	Infants 5-10 kg	M0351	1 each
2	Infants/Children 10-20 kg	M0361	1 each
2½	Children 20-30 kg	M0371	1 each
3	Children 30-50 kg	M0311	1 each
4	Adults 50-70 kg	M0321	1 each
5	Adults 70-100 kg	M0331	1 each
6	Adults > 100 kg	M0381	1 each

PEDIATRIC



SourceFlex™ Reusable Laryngeal Mask Airway



FEATURES AND BENEFITS

- Reinforced with wire to resist kinking and maintain cuff seal during manipulation of patient
- Used in ENT, head & neck and upper body procedures
- Multi-use convenience
- 100% medical grade silicone with ABS/PC connector/valve
- No aperture bars for easier insertion of diagnostic tools and cleaning

MASK SIZE	DESCRIPTION	ITEM #	QTY	PEDIATRIC
2	Infants/Children 10-20 kg	M0364	1 each	
2½	Children 20-30 kg	M0374	1 each	
3	Children 30-50 kg	M0314	1 each	
4	Adults 50-70 kg	M0324	1 each	
5	Adults 70-100 kg	M0334	1 each	

SourceFlex™ Disposable Silicone Laryngeal Mask Airway



FEATURES AND BENEFITS

- Reinforced with wire to resist kinking and maintain cuff seal during manipulation of patient
- 100% medical grade silicone with ABS/PC connector/valve
- Color coded pilot balloon for quick size identification
- No aperture bars for easier insertion of diagnostic tools
- Used in ENT, head & neck and upper body procedures

MASK SIZE	DESCRIPTION	ITEM #	QTY	PEDIATRIC
2	Infants/Children 10-20 kg	M0365	5 / box	
2½	Children 20-30 kg	M0375	5 / box	
3	Children 30-50 kg	M0315	5 / box	
4	Adults 50-70 kg	M0325	5 / box	
5	Adults 70-100 kg	M0335	5 / box	



Laryngeal Mask Airway with Gastric Access






FEATURES AND BENEFITS

- Integrated gastric access channel
- Anatomical curve for easy insertion
- Built in bite block to prevent compromised airway
- Color coded pilot balloon for quick size identification
- Can be used with OG tubes up to size 14FR

Intended Use

- Abdominal procedures
- Patients at risk of massive reflux
- Unexpected difficult airways
- Patients who have not fasted

MASK SIZE	DESCRIPTION	ITEM #	QTY
 3	Children 30-50 kg	M0316	5 / box
 4	Adults 50-70 kg	M0326	5 / box
 5	Adults 70-100 kg	M0336	5 / box



Enhanced Comfort Endotracheal Tube



FEATURES AND BENEFITS

- DEHP Free
- High volume, low profile cylindrical cuff
- Made of softer PVC material, reducing risk of trauma

TUBE SIZE	ITEM #	QTY
5mm	M0450CS	10/box
5.5mm	M0455CS	10/box
6mm	M0460CS	10/box
6.5mm	M0465CS	10/box
7mm	M0470CS	10/box
7.5mm	M0475CS	10/box
8mm	M0480CS	10/box
8.5mm	M0485CS	10/box
9mm	M0490CS	10/box

SuctionPlus ® Endotracheal Tube



FEATURES AND BENEFITS

- Latex & DEHP Free
- Removable standard 15mm connector
- Thermosmoothing of the tube tips
- High volume, low pressure cylindrical cuff
- One-way inflation valve with inflation pilot balloon
- Sizes 6.0-9.0mm, individual sterile package

Intended Use

SuctionPlus Endotracheal tube is designed with the evacuation lumen in order to provide protection against risk of aspiration and decrease the rate of ventilation-associated pneumonia (VAP) in the hospital.


TUBE SIZE	ITEM #	QTY
6mm	M0495C	10/box
6.5mm	M0500C	10/box
7mm	M0505C	10/box
7.5mm	M0510C	10/box
8mm	M0515C	10/box
8.5mm	M0520C	10/box
9mm	M0525C	10/box



Standard Hi-Lo Cuffed ET Tube

10 ea per box / 10 boxes per case

Hi/Lo Straight HVLP Endotracheal Tube, Cuffed, Murphy, 15mm connector



Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
ITEM #	M0430C	M0435C	M0440C	M0445C	M0450C	M0455C	M0460C	M0465C	M0470C	M0475C	M0480C	M0485C	M0490C

Oral Cuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Oral Preformed, Cuffed, Murphy, 15mm connector



Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5
ITEM #	M0530C	M0535C	M0540C	M0545C	M0550C	M0555C	M0560C	M0565C	M0570C	M0575C	M0580C	M0585C

Nasal Cuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Nasal Preformed, Cuffed, Murphy, 15mm connector



Tube Size	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0
ITEM #	M0640C	M0645C	M0650C	M0655C	M0660C	M0665C	M0670C	M0675C	M0680C

Reinforced Cuffed ET Tube

5 ea per box / 10 boxes per case

Endotracheal Tube, Reinforced, Cuffed, Murphy, 15mm connector



Tube Size	5.5	6.0	6.5	7.0	7.5	8.0	8.5
ITEM #	M0755C	M0760C	M0765C	M0770C	M0775C	M0780C	M0785C

The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

CUFFED:

- High-volume low-pressure cylindrical cuff distributes even pressure on the tracheal wall reducing chance of tissue damage
- Large cuff surface provides better seal
- Durable construction to ensure the reliability of the tube and cuff during procedures

TUBE:

- Thermo-smoothing of the tube tips to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized



Standard Uncuffed ET Tube

10 ea per box / 10 boxes per case

Straight Endotracheal Tube, Uncuffed, Murphy, 15mm connector

Tube Size	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0
ITEM #	M0420U	M0425U	M0430U	M0435U	M0440U	M0445U	M0450U	M0455U	M0460U



Oral Uncuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Oral Preformed, Uncuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0
ITEM #	M0530U	M0535U	M0540U	M0545U	M0550U	M0555U	M0560U



Nasal Uncuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Nasal Preformed, Uncuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0
ITEM #	M0630U	M0635U	M0640U	M0645U	M0650U	M0655U	M0660U	M0665U	M0670U



Reinforced Uncuffed ET Tube

5 ea per box / 10 boxes per case

Endotracheal Tube, Reinforced, Cuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0
ITEM #	M0730U	M0735U	M0740U	M0745U	M0750U



The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

TUBE:

- Thermo-smoothing of the tube tips to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized




Standard Hi-Lo Cuffed ET Tube with Pre-Loaded Stylet

10 ea per box / 10 boxes per case

Tubes, Endo, Oral and Nasal, Murphy, Cuffed + Stylet, Aluminum

Tube Size	6.5	7.0	7.5	8.0
ITEM #	MX0465C	MX0470C	MX0475C	MX0480C

SourceMark's Hi-Lo cuffed Endotracheal Tube with Stylet preloaded for quick intubation



Aluminum Stylet

20 ea per box / 20 boxes per case

Indications for usage:

- Airway management during procedures
- Anesthesia

FEATURES AND BENEFITS

- Malleable aluminum which allows the stylet to be easily adapted to any desired shape
- Coated to reduce friction between the stylet and the tracheal tube to allow for easy insertion and withdrawal

Size	6FR	8FR	10FR	12FR	14FR
ITEM #	M0906	M0908	M0910	M0912	M0914



Nasopharyngeal Airways

10 ea per box / 10 boxes per case

FEATURES AND BENEFITS

- Maximum kink resistance
- Smoothly rounded edges for patient comfort
- Made of medical grade PVC
- Angled design to reduce risk of injury
- Beveled tip for easy insertion
- Available in full range of sizes

Size	14FR	16FR	18FR	20FR	22FR	24FR	26FR	28FR	30FR	32FR	34FR	36FR
ITEM #	M1014	M1016	M1018	M1020	M1022	M1024	M1026	M1028	M1030	M1032	M1034	M1036





ClearViz™ Single-Use Conventional Laryngoscope Blades

FEATURES AND BENEFITS

- Enclosed fiber optic cable
- Zinc metal blade construction
- Matte Finish
- Compatible with ClearViz orange handle

Item #	Description	Qty
SLMAC10	Macintosh Blade Size 1	10/box
SLMAC20	Macintosh Blade Size 2	10/box
SLMAC30	Macintosh Blade Size 3	10/box
SLMAC40	Macintosh Blade Size 4	10/box



Item #	Description	Qty
SLMIL00	Miller Blade Size 00	10/box
SLMIL0	Miller Blade Size 0	10/box
SLMIL10	Miller Blade Size 1	10/box
SLMIL20	Miller Blade Size 2	10/box
SLMIL30	Miller Blade Size 3	10/box
SLMIL40	Miller Blade Size 4	10/box



ClearViz™ Single-Use Conventional Laryngoscope Handle

FEATURES AND BENEFITS

- Textured Grip Handle
- Lightweight and Durable
- Latex & DEHP free
- LED 1W bulb with fiberoptic cable generates >2K lumens

Item #	Description	Qty
SLH02	Laryngoscope Handle, Disposable, Medium	10/box





Esophageal Stethoscope

FEATURES AND BENEFITS

- Dual purpose device combining internal auscultation of heart and breath sounds with temperature sensor for monitoring patient
- Placed in the lower esophagus at the heart level providing an auditory transmission of precise heart and breath sounds directly to the anesthesiologist's earpiece
- Specially formulated balloon material to optimize acoustics and sound transmission
- Balloon heat bonded to tube with smooth finish to minimize risk of esophageal trauma during insertion and removal
- Long lead wires keep connections off patient's face and out of the way of the anesthesiologist during surgery
- Available in multiple sizes for use in infant to adult

Item #	Description	Qty
SM9400ES	Esophageal Stethoscope w/ 400 Series Thermistor 9F	30/box
SM12400ES	Esophageal Stethoscope w/ 400 Series Thermistor 12F	30/box
SM18400ES	Esophageal Stethoscope w/ 400 Series Thermistor 18F	30/box
SM24400ES	Esophageal Stethoscope w/ 400 Series Thermistor 24F	30/box



General Purpose Probe

FEATURES AND BENEFITS

- Provides accurate measurement of body temperature when placed in the esophagus, rectum, nasopharynx or axilla
- Tip-to-end PVC sleeve protects the lead wire from strain, breaks, kinks and interference from ambient air, enhancing product quality and temperature accuracy
- Smooth, rounded beveled tip for non-traumatic insertion and withdrawal to the patient
- Available in 9FR and 12FR sizes for use in infant to adult
- Compatible with all Series 400 patient monitoring equipment and all hypo-hyperthermia equipment providing an accurate, affordable alternative

Item #	Description	Qty
SM9400ER	General Purpose Esophageal/Rectal Probe w/ 400 Series Thermistor 9F	35/box
SM12400ER	General Purpose Esophageal/Rectal Probe w/ 400 Series Thermistor 12F	35/box



Skin Temperature Sensor

FEATURES AND BENEFITS

- Sensor combines a precision thermistor with soft, pliable, adhesive foam disk for easy and secure attachment directly to patient's skin
- Sensor is securely embedded into the center of the foam disk to reduce pressure and irritation on the skin and minimize any risk of sensor separation from the foam
- Thinner lead wire offers more pliable application on the skin and reduces the interference of ambient air on the temperature reading (the stem effect)
- Mylar backing provides insulation from heat lamps and other lighting, improving sensor accuracy

Item #	Description	Qty
SM400STS	Skin Temperature Sensor w/ 400 Series Thermistor	50/box





Neonatal Temperature Probe

FEATURES AND BENEFITS

- Thermal Reflectors hold Temperature Probes securely to skin
- Hydrogel is gentle on infant skin and allows easy repositioning
- Temperature Probes are single-patient use for length of stay
- Temperature Probes come with 100 Large Reflectors

Item #	Description	Qty
SM4200	Thermal Reflector, 1.5in X 1.25in, Large, Gold Heart-Shaped	200/case
SM4220	Thermal Reflector, 1in X 7/8in, Small, Gold Heart-Shaped	100/case
SM4501NP	Neonatal Temp Probes for Air-Shields® Isolette® Infant Incubator Model C-100. Infant Radiant Warmer 7810, 7820, 7830, 7840, 7860, 7865	50/case
SM4503NP	Neonatal Temp Probes for Ohmeda® IC, GC, and CP Incubators. Infant Warmer Systems 3000, 3300, and 5000	50/case
SM4504NP	Neonatal Temp Probes for GE® Giraffe OmniBed, Giraffe Incubators, Giraffe Warmers	50/case
SM4505NP	Neonatal Temp Probes for Dräger® Babytherm® 8004 & Babytherm®. 8010 CALEO Incubators	50/case
SM4506NP	Neonatal Temp Probes for Dräger® Isolette® C2000 & 8000 Incubators. Versalet 7700 Care Center	50/case
SM4507NP	Neonatal Temp Probes for ATOM Infa Warmer 1, ATOM Medical Incubators	50/case



Reflectors



SM4501NP



SM4503NP



SM4504NP



SM4505NP



SM4506NP



SM4507NP

Tympanic Temperature Probe

FEATURES AND BENEFITS

- Provides safe and secure seal of ear canal for canal to equilibrate to tympanic temperature
- Probe features a dual memory foam tip for protection of the tympanic membrane with less risk of perforation or damage

Item #	Description	Qty
SM400TY	Adult Tympanic Temperature Sensor 400 Series	50/box
SM400TYP	Pediatric Tympanic Temperature Sensor 400 Series	50/box



Temperature Probe Cable

Item #	Description	Qty
SMC400-MP-M	Standard 1/4" Phone plug 400 Series, 10ft	10/box
SMC400-P-M	Philips/HP Cable 400 Series, 10ft	10/box





SourceMarkTM
Medical