



Tourniquet Cuffs

Phone: 615.269.6010
Toll Free: 866.528.7001
Fax: 866.528.6008

www.SourceMarkUSA.com

TOURNIQUET CUFFS

Tourniquet Dispo Cuff / Single Cuff / Single Port

Tourniquet cuff for single use - sterile



Length	Width	Colour	Single Port	Box
20 cm (8")	5 cm (2")	-	REF 20-34-700SLZ-1	10
30 cm (12")	6 cm (2 3/8")	Blue	REF 20-34-710SLZ-1	10
35 cm (14")	8 cm (3 1/8")	Yellow	REF 20-34-711SLZ-1	10
46 cm (18")	10 cm (4")	Red	REF 20-34-712SLZ-1	10
46 cm (18")	12 cm (4 3/4")	Red	REF 20-34-715SLZ-1*	10
61 cm (24")	12 cm (4 3/4")	Green	REF 20-34-722SLZ-1*	10
76 cm (30")	12 cm (4 3/4")	Dark Blue	REF 20-34-727SLZ-1*	10
86 cm (34")	12 cm (4 3/4")	Brown	REF 20-34-728SLZ-1*	10
107 cm (42")	12 cm (4 3/4")	Dark Red	REF 20-34-729SLZ-1*	10

*Contour shape



Tourniquet Dispo Cuff / Single Cuff / Dual Port

Tourniquet cuff for single use - sterile



Length	Width	Colour	Dual Port	Box
20 cm (8")	5 cm (2")	-	REF 20-36-700SLZ-1	10
30 cm (12")	6 cm (2 3/8")	Blue	REF 20-36-710SLZ-1	10
35 cm (14")	8 cm (3 1/8")	Yellow	REF 20-36-711SLZ-1	10
46 cm (18")	10 cm (4")	Red	REF 20-36-712SLZ-1	10
46 cm (18")	12 cm (4 3/4")	Red	REF 20-36-715SLZ-1*	10
61 cm (24")	12 cm (4 3/4")	Green	REF 20-36-722SLZ-1*	10
76 cm (30")	12 cm (4 3/4")	Dark Blue	REF 20-36-727SLZ-1*	10
86 cm (34")	12 cm (4 3/4")	Brown	REF 20-36-728SLZ-1*	10
107 cm (42")	12 cm (4 3/4")	Dark Red	REF 20-36-729SLZ-1*	10

*Contour shape

TOURNIQUET CUFFS

Tourniquet Dispo Cuff / Double Cuff / Single Port

Tourniquet cuff for single use - sterile



Length	Width	Colour	REF	Box
30 cm (12")	12 cm (4 3/4")	Blue	20-30-710SLZ-1	10
46 cm (18")	12 cm (6")	Red	20-30-712SLZ-1	10
61 cm (24")	15 cm (6")	Green	20-30-722SLZ-1	10

Tourniquet Wipe Cuff / Single Cuff / Single Port

Tourniquet cuff - reusable




Length	Width	Colour	REF	Box
20 cm (8")	7 cm (2 3/4")	-	20-75-700	1
30 cm (12")	8 cm (3 1/8")	Blue	20-75-710	1
35 cm (14")	8 cm (3 1/8")	Yellow	20-75-711	1
46 cm (18")	10 cm (4")	Red	20-75-712	1
46 cm (18")	12 cm (4 3/4")	Red	20-75-715*	1
61 cm (24")	12 cm (4 3/4")	Green	20-75-722*	1
76 cm (30")	12 cm (4 3/4")	Dark Blue	20-75-727*	1
86 cm (34")	12 cm (4 3/4")	Brown	20-75-728*	1
107 cm (42")	12 cm (4 3/4")	Dark Red	20-75-729*	1

*Contour shape

Tourniquet Wipe Cuff / Double Cuff / Single Port

Tourniquet cuff - reusable



Length	Width	Colour	REF	Box
30 cm (12")	12 cm (4.7")	Blue	20-77-710	1
46 cm (18")	15 cm (5.9")	Red	20-77-712	1
61 cm (24")	15 cm (5.9")	Green	20-77-722	1

TOURNIQUET CUFFS


Silicone Velcro Cuff / Single Cuff / Single Port

Tourniquet cuff - reusable



Length	Width	Colour	REF	Box
20 cm (8")	2.5 cm (1")	-	20-64-700	1
30 cm (12")	4 cm (1 5/8")	Blue	20-64-710	1
35 cm (14")	6 cm (2 3/8")	Yellow	20-64-611	1
46 cm (18")	7.5 cm (3")	Red	20-64-612	1
46 cm (18")	11 cm (4 1/4")	Red	20-64-512*	1
61 cm (24")	11 cm (4 1/4")	Green	20-64-522*	1
76 cm (30")	11 cm (4 1/4")	Dark Blue	20-64-527*	1
86 cm (34")	11 cm (4 1/4")	Brown	20-64-528*	1

*Contour shape



Silicone Velcro Cuff / Double Cuff / Single Port

Tourniquet cuff - reusable



Length	Width	Colour	REF	Box
30 cm (12")	8.0 cm (3.2")	Blue	20-60-710	1
46 cm (18")	12.5 cm (4.9")	Red	20-60-712	1
61 cm (24")	12.5 cm (4.9")	Green	20-60-722	1

TOURNIQUET CUFFS

Disposable Tourniquet Cuff with Sleeve / Single Cuff / Single Port

Tourniquet cuff for single use - sterile



Part Number	Size	Colour	
M8806SP	46 cm (18")	■	Standard Arm, Single Port, with Sleeve
M8808SP	61 cm (24")	■	Standard Leg, Single Port, with Sleeve
M8810SP	76 cm (30")	□	Medium Leg, Single Port, with Sleeve
M8812SP	86 cm (34")	■	Long Leg, Single Port, with Sleeve
M8814SP	111 cm (44")	■	Extra Long Leg, Single Port, with Sleeve

Disposable Tourniquet Cuff with Sleeve / Single Cuff / Dual Port

Tourniquet cuff for single use - sterile



Part Number	Size	Colour	
M8807DP	46 cm (18")	■	Standard Arm, Dual Port, with Sleeve
M8809DP	61 cm (24")	■	Standard Leg, Dual Port, with Sleeve
M8811DP	76 cm (30")	□	Medium Leg, Dual Port, with Sleeve
M8813DP	86 cm (34")	■	Long Leg, Dual Port, with Sleeve
M8815DP	111 cm (44")	■	Extra Long Leg, Dual Port, with Sleeve





SourceMark™
Medical