



Intubating Video Laryngeal Mask Airway

Seeing is believing ...

Surveillance of the supraglottic airway



Benefits:

1. Precise insertion provides rapid and reliable ventilation
2. Intraoperative monitoring ensures airway safety
3. Video-assisted tracheal intubation and extubation


SOURCEMARK MEDICAL™

📍 302 INNOVATION DRIVE, SUITE 410, FRANKLIN TN, 37067

✉ CUSTOMERSERVICE@SOURCEMARKUSA.COM

🌐 WWW.SOURCEMARKUSA.COM

Supraglottic airway with a vision™



Product Specification

Model	Size	Patient Weight (kg)	Gastric Tube (Fr)	Max. ETT Size (mm)	Video Stylet Model
LMUE010M	1.0	< 5	6	3.5	VS-10
LMUE015M	1.5	5-10	8	4.0	VS-10
LMUE020M	2.0	10-20	10	5.0	VS-20
LMUE025M	2.5	20-30	10	5.5	VS-20
LMUE030M	3.0	30-50	16	7.0	VS-40
LMUE040M	4.0	50-70	16	7.5	VS-40
LMUE050M	5.0	70-100	16	8.0	VS-40