



Available as individual blades and handles, as well as convenient blade and handle combo packs!

Part Number	Combo Packs
SLMCH10	Mac Blade & Handle, Disposable, Size 1.0
SLMCH20	Mac Blade & Handle, Disposable, Size 2.0
SLMCH30	Mac Blade & Handle, Disposable, Size 3.0
SLMCH40	Mac Blade & Handle, Disposable, Size 4.0
SLMLH0	Miller Blade & Handle, Disposable, Size 0
SLMLH00	Miller Blade & Handle, Disposable, Size 00
SLMLH10	Miller Blade & Handle, Disposable, Size 1.0
SLMLH20	Miller Blade & Handle, Disposable, Size 2.0
SLMLH30	Miller Blade & Handle, Disposable, Size 3.0
SLMLH40	Miller Blade & Handle, Disposable, Size 4.0



ClearViz™

Blade and Handle Details

ClearViz Blade Advantage

- Durable zinc metal blade construction
- Fully encapsulated fiber optics to reduce light refraction
- Large easy to read blade size marking for quick identification
- Matte coat blade finish to minimize light reflection
- Color coded compatibility for rapid pairing

ClearViz Handle Advantage

- Robust and lightweight
- Textured grip handle
- Super bright LED- >2K lumens
- Powered by 2ea AAA batteries
- Latex & DEHP free

Part Number	Individual Blades
SLMAC10	Macintosh Blade Size 1
SLMAC20	Macintosh Blade Size 2
SLMAC30	Macintosh Blade Size 3
SLMAC40	Macintosh Blade Size 4



Part Number	Individual Blades
SLMIL00	Miller Blade Size 00
SLMILO	Miller Blade Size 0
SLMIL10	Miller Blade Size 1
SLMIL20	Miller Blade Size 2
SLMIL30	Miller Blade Size 3
SLMIL40	Miller Blade Size 4



Part Number	Individual Handle
SLH02	Laryngoscope Handle, Disposable, Medium

