



SourceMark™
Medical



ANESTHESIA / AIRWAY MANAGEMENT



Intubating Video Laryngeal Mask Airway



MASK SIZE	DESCRIPTION	ITEM #	VIDEO STYLET MODEL
1	Neonates/Infants up to 5 kg	MLM10C	VS-10
1½	Infants 5-10 kg	MLM15C	VS-10
2	Infants/Children 10-20 kg	MLM20C	VS-20
2½	Children 20-30 kg	MLM25C	VS-20
3	Children 30-50 kg	MLM30C	VS-40
4	Adults 50-70 kg	MLM40C	VS-40
5	Adults 70-100 kg	MLM50C	VS-40

PEDIATRIC

FEATURES AND BENEFITS

- Precise insertion provides rapid and reliable ventilation
- Intraoperative monitoring ensures airway safety
- Enables video-assisted tracheal intubation and extubation
- Rapid ventilation for normal and difficult airways
- The preferred choice for rescue ventilation in various clinical scenarios

Supraglottic airway with a vision

See-through camera port with cleaning technology

End of gastric channel

Inflation line

Inflatable cuff

Cleaning port to flush camera lens


Observation channel
Ventilation (intubation) channel
Gastric channel
Bite block

Demo Video





Silicone Disposable Laryngeal Mask Airway



FEATURES AND BENEFITS

- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, neonate through adult
- Color coded pilot balloon for quick size identification
- No aperture bars for easier insertion of diagnostic tools
- Supple cuff material for less irritation

MASK SIZE	DESCRIPTION	ITEM #	QTY
1	Neonates/Infants up to 5 kg	M0342	10/box
1½	Infants 5-10 kg	M0352	10/box
2	Infants/Children 10-20 kg	M0362	10/box
2½	Children 20-30 kg	M0372	10/box
3	Children 30-50 kg	M0312	10/box
4	Adults 50-70 kg	M0322	10/box
5	Adults 70-100 kg	M0332	10/box

PEDIATRIC

Standard Disposable Laryngeal Mask Airway



FEATURES AND BENEFITS

- No aperture bars for easier insertion of diagnostic tools
- Anatomically correct curve
- Sturdy tube for easier insertion
- Easy product transition
- Available in full range of sizes, neonate through adult

MASK SIZE	DESCRIPTION	ITEM #	QTY
1	Neonates/Infants up to 5 kg	M0343	10/box
1½	Infants 5-10 kg	M0353	10/box
2	Infants/Children 10-20 kg	M0363	10/box
2½	Children 20-30 kg	M0373	10/box
3	Children 30-50 kg	M0313	10/box
4	Adults 50-70 kg	M0323	10/box
5	Adults 70-100 kg	M0333	10/box

PEDIATRIC



Quikpak Silicone Disposable Laryngeal Mask Airway



FEATURES AND BENEFITS

- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, infants through adult
- Color coded pilot balloon for quick size identification
- No aperture bars for easier insertion of diagnostic tools
- Supple cuff material for less irritation

MASK SIZE	DESCRIPTION	ITEM #	QTY
1½	Infants 5-10 kg	MX0352	10/box
2	Infants/Children 10-20 kg	MX0362	10/box
2½	Children 20-30 kg	MX0372	10/box
3	Children 30-50 kg	MX0312	10/box
4	Adults 50-70 kg	MX0322	10/box
5	Adults 70-100 kg	MX0332	10/box

PEDIATRIC

*QuikPak is SourceMark's MX series for fast and convenient usage.
It includes one syringe and one pack of 5.0g (.17oz) of lubricating jelly.

Reusable Laryngeal Mask Airway



FEATURES AND BENEFITS


- 100% medical grade silicone with ABS/PC connector/valve
- Available in full range of sizes, neonate through adult
- Multi-use convenience
- No aperture bars for easier insertion of diagnostic tools and cleaning
- Supple cuff material for less irritation

MASK SIZE	DESCRIPTION	ITEM #	QTY
1	Neonates/Infants up to 5 kg	M0341	1 each
1½	Infants 5-10 kg	M0351	1 each
2	Infants/Children 10-20 kg	M0361	1 each
2½	Children 20-30 kg	M0371	1 each
3	Children 30-50 kg	M0311	1 each
4	Adults 50-70 kg	M0321	1 each
5	Adults 70-100 kg	M0331	1 each
6	Adults > 100 kg	M0381	1 each

PEDIATRIC



SourceFlex™ Reusable Laryngeal Mask Airway



FEATURES AND BENEFITS

- Reinforced with wire to resist kinking and maintain cuff seal during manipulation of patient
- Used in ENT, head & neck and upper body procedures
- Multi-use convenience
- 100% medical grade silicone with ABS/PC connector/valve
- No aperture bars for easier insertion of diagnostic tools and cleaning

MASK SIZE	DESCRIPTION	ITEM #	QTY	PEDIATRIC
2	Infants/Children 10-20 kg	M0364	1 each	
2½	Children 20-30 kg	M0374	1 each	
3	Children 30-50 kg	M0314	1 each	
4	Adults 50-70 kg	M0324	1 each	
5	Adults 70-100 kg	M0334	1 each	

SourceFlex™ Disposable Silicone Laryngeal Mask Airway



FEATURES AND BENEFITS

- Reinforced with wire to resist kinking and maintain cuff seal during manipulation of patient
- 100% medical grade silicone with ABS/PC connector/valve
- Color coded pilot balloon for quick size identification
- No aperture bars for easier insertion of diagnostic tools
- Used in ENT, head & neck and upper body procedures

MASK SIZE	DESCRIPTION	ITEM #	QTY	PEDIATRIC
2	Infants/Children 10-20 kg	M0365	5 / box	
2½	Children 20-30 kg	M0375	5 / box	
3	Children 30-50 kg	M0315	5 / box	
4	Adults 50-70 kg	M0325	5 / box	
5	Adults 70-100 kg	M0335	5 / box	



Laryngeal Mask Airway with Gastric Access






FEATURES AND BENEFITS

- Integrated gastric access channel
- Anatomical curve for easy insertion
- Built in bite block to prevent compromised airway
- Color coded pilot balloon for quick size identification
- Can be used with OG tubes up to size 14FR

Intended Use

- Abdominal procedures
- Patients at risk of massive reflux
- Unexpected difficult airways
- Patients who have not fasted

MASK SIZE	DESCRIPTION	ITEM #	QTY
 3	Children 30-50 kg	M0316	5 / box
 4	Adults 50-70 kg	M0326	5 / box
 5	Adults 70-100 kg	M0336	5 / box



Enhanced Comfort Endotracheal Tube



FEATURES AND BENEFITS

- DEHP Free
- High volume, low profile cylindrical cuff
- Made of softer PVC material, reducing risk of trauma

TUBE SIZE	ITEM #	QTY
5mm	M0450CS	10/box
5.5mm	M0455CS	10/box
6mm	M0460CS	10/box
6.5mm	M0465CS	10/box
7mm	M0470CS	10/box
7.5mm	M0475CS	10/box
8mm	M0480CS	10/box
8.5mm	M0485CS	10/box
9mm	M0490CS	10/box

SuctionPlus ® Endotracheal Tube



FEATURES AND BENEFITS

- Latex & DEHP Free
- Removable standard 15mm connector
- Thermosmoothing of the tube tips
- High volume, low pressure cylindrical cuff
- One-way inflation valve with inflation pilot balloon
- Sizes 6.0-9.0mm, individual sterile package

Intended Use

SuctionPlus Endotracheal tube is designed with the evacuation lumen in order to provide protection against risk of aspiration and decrease the rate of ventilation-associated pneumonia (VAP) in the hospital.


TUBE SIZE	ITEM #	QTY
6mm	M0495C	10/box
6.5mm	M0500C	10/box
7mm	M0505C	10/box
7.5mm	M0510C	10/box
8mm	M0515C	10/box
8.5mm	M0520C	10/box
9mm	M0525C	10/box



Standard Hi-Lo Cuffed ET Tube

10 ea per box / 10 boxes per case

Hi/Lo Straight HVLP Endotracheal Tube, Cuffed, Murphy, 15mm connector



Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
ITEM #	M0430C	M0435C	M0440C	M0445C	M0450C	M0455C	M0460C	M0465C	M0470C	M0475C	M0480C	M0485C	M0490C

Oral Cuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Oral Preformed, Cuffed, Murphy, 15mm connector



Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5
ITEM #	M0530C	M0535C	M0540C	M0545C	M0550C	M0555C	M0560C	M0565C	M0570C	M0575C	M0580C	M0585C

Nasal Cuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Nasal Preformed, Cuffed, Murphy, 15mm connector



Tube Size	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0
ITEM #	M0640C	M0645C	M0650C	M0655C	M0660C	M0665C	M0670C	M0675C	M0680C

Reinforced Cuffed ET Tube

5 ea per box / 10 boxes per case

Endotracheal Tube, Reinforced, Cuffed, Murphy, 15mm connector



Tube Size	5.5	6.0	6.5	7.0	7.5	8.0	8.5
ITEM #	M0755C	M0760C	M0765C	M0770C	M0775C	M0780C	M0785C

The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

CUFFED:

- High-volume low-pressure cylindrical cuff distributes even pressure on the tracheal wall reducing chance of tissue damage
- Large cuff surface provides better seal
- Durable construction to ensure the reliability of the tube and cuff during procedures

TUBE:

- Thermo-smoothing of the tube tips to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized



Standard Uncuffed ET Tube

10 ea per box / 10 boxes per case

Straight Endotracheal Tube, Uncuffed, Murphy, 15mm connector

Tube Size	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0
ITEM #	M0420U	M0425U	M0430U	M0435U	M0440U	M0445U	M0450U	M0455U	M0460U



Oral Uncuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Oral Preformed, Uncuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0
ITEM #	M0530U	M0535U	M0540U	M0545U	M0550U	M0555U	M0560U



Nasal Uncuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Nasal Preformed, Uncuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0
ITEM #	M0630U	M0635U	M0640U	M0645U	M0650U	M0655U	M0660U	M0665U	M0670U



Reinforced Uncuffed ET Tube

5 ea per box / 10 boxes per case

Endotracheal Tube, Reinforced, Cuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0
ITEM #	M0730U	M0735U	M0740U	M0745U	M0750U



The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

TUBE:

- Thermo-smoothing of the tube tips to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized




Standard Hi-Lo Cuffed ET Tube with Pre-Loaded Stylet

10 ea per box / 10 boxes per case

Tubes, Endo, Oral and Nasal, Murphy, Cuffed + Stylet, Aluminum

Tube Size	6.5	7.0	7.5	8.0
ITEM #	MX0465C	MX0470C	MX0475C	MX0480C

SourceMark's Hi-Lo cuffed Endotracheal Tube with Stylet preloaded for quick intubation



Aluminum Stylet

20 ea per box / 20 boxes per case

Indications for usage:

- Airway management during procedures
- Anesthesia

FEATURES AND BENEFITS

- Malleable aluminum which allows the stylet to be easily adapted to any desired shape
- Coated to reduce friction between the stylet and the tracheal tube to allow for easy insertion and withdrawal

Size	6FR	8FR	10FR	12FR	14FR
ITEM #	M0906	M0908	M0910	M0912	M0914



Nasopharyngeal Airways

10 ea per box / 10 boxes per case

FEATURES AND BENEFITS

- Maximum kink resistance
- Smoothly rounded edges for patient comfort
- Made of medical grade PVC
- Angled design to reduce risk of injury
- Beveled tip for easy insertion
- Available in full range of sizes

Size	14FR	16FR	18FR	20FR	22FR	24FR	26FR	28FR	30FR	32FR	34FR	36FR
ITEM #	M1014	M1016	M1018	M1020	M1022	M1024	M1026	M1028	M1030	M1032	M1034	M1036





ClearViz™ Single-Use Conventional Laryngoscope Blades and Handles

FEATURES AND BENEFITS

- Enclosed fiber optic cable
- Zinc metal blade construction
- Matte Finish
- Compatible with ClearViz orange handle

Item #	Individual Blades	Qty
SLMAC10	Macintosh Blade Size 1	10/box
SLMAC20	Macintosh Blade Size 2	10/box
SLMAC30	Macintosh Blade Size 3	10/box
SLMAC40	Macintosh Blade Size 4	10/box



Item #	Individual Blades	Qty
SLMIL00	Miller Blade Size 00	10/box
SLMIL0	Miller Blade Size 0	10/box
SLMIL10	Miller Blade Size 1	10/box
SLMIL20	Miller Blade Size 2	10/box
SLMIL30	Miller Blade Size 3	10/box
SLMIL40	Miller Blade Size 4	10/box



ClearViz™ Single-Use Conventional Laryngoscope Handle

FEATURES AND BENEFITS

- Textured Grip Handle
- Lightweight and Durable
- Latex & DEHP free
- LED 1W bulb with fiberoptic cable generates >2K lumens

Item #	Individual Handles	Qty
SLH02	Laryngoscope Handle, Disposable, Medium	10/box





ClearViz™ Single-Use Conventional Laryngoscope Combo Packs

Item #	Combo Packs	Qty
SLMCH10	Mac Blade & Handle, Disposable, Size 1.0	8/box
SLMCH20	Mac Blade & Handle, Disposable, Size 2.0	8/box
SLMCH30	Mac Blade & Handle, Disposable, Size 3.0	8/box
SLMCH40	Mac Blade & Handle, Disposable, Size 4.0	8/box
SLMLH0	Miller Blade & Handle, Disposable, Size 0	8/box
SLMLH00	Miller Blade & Handle, Disposable, Size 00	8/box
SLMLH10	Miller Blade & Handle, Disposable, Size 1.0	8/box
SLMLH20	Miller Blade & Handle, Disposable, Size 2.0	8/box
SLMLH30	Miller Blade & Handle, Disposable, Size 3.0	8/box
SLMLH40	Miller Blade & Handle, Disposable, Size 4.0	8/box





Esophageal Stethoscope

FEATURES AND BENEFITS

- Dual purpose device combining internal auscultation of heart and breath sounds with temperature sensor for monitoring patient
- Placed in the lower esophagus at the heart level providing an auditory transmission of precise heart and breath sounds directly to the anesthesiologist's earpiece
- Specially formulated balloon material to optimize acoustics and sound transmission
- Balloon heat bonded to tube with smooth finish to minimize risk of esophageal trauma during insertion and removal
- Long lead wires keep connections off patient's face and out of the way of the anesthesiologist during surgery
- Available in multiple sizes for use in infant to adult

Item #	Description	Qty
SM9400ES	Esophageal Stethoscope w/ 400 Series Thermistor 9F	30/box
SM12400ES	Esophageal Stethoscope w/ 400 Series Thermistor 12F	30/box
SM18400ES	Esophageal Stethoscope w/ 400 Series Thermistor 18F	30/box
SM24400ES	Esophageal Stethoscope w/ 400 Series Thermistor 24F	30/box



General Purpose Probe

FEATURES AND BENEFITS

- Provides accurate measurement of body temperature when placed in the esophagus, rectum, nasopharynx or axilla
- Tip-to-end PVC sleeve protects the lead wire from strain, breaks, kinks and interference from ambient air, enhancing product quality and temperature accuracy
- Smooth, rounded beveled tip for non-traumatic insertion and withdrawal to the patient
- Available in 9FR and 12FR sizes for use in infant to adult
- Compatible with all Series 400 patient monitoring equipment and all hypo-hyperthermia equipment providing an accurate, affordable alternative

Item #	Description	Qty
SM9400ER	General Purpose Esophageal/Rectal Probe w/ 400 Series Thermistor 9F	35/box
SM12400ER	General Purpose Esophageal/Rectal Probe w/ 400 Series Thermistor 12F	35/box



Skin Temperature Sensor

FEATURES AND BENEFITS

- Sensor combines a precision thermistor with soft, pliable, adhesive foam disk for easy and secure attachment directly to patient's skin
- Sensor is securely embedded into the center of the foam disk to reduce pressure and irritation on the skin and minimize any risk of sensor separation from the foam
- Thinner lead wire offers more pliable application on the skin and reduces the interference of ambient air on the temperature reading (the stem effect)
- Mylar backing provides insulation from heat lamps and other lighting, improving sensor accuracy

Item #	Description	Qty
SM400STS	Skin Temperature Sensor w/ 400 Series Thermistor	50/box





Neonatal Temperature Probe

FEATURES AND BENEFITS

- Thermal Reflectors hold Temperature Probes securely to skin
- Hydrogel is gentle on infant skin and allows easy repositioning
- Temperature Probes are single-patient use for length of stay
- Temperature Probes come with 100 Large Reflectors

Item #	Description	Qty
SM4200	Thermal Reflector, 1.5in X 1.25in, Large, Gold Heart-Shaped	200/case
SM4220	Thermal Reflector, 1in X 7/8in, Small, Gold Heart-Shaped	100/case
SM4501NP	Neonatal Temp Probes for Air-Shields® Isolette® Infant Incubator Model C-100. Infant Radiant Warmer 7810, 7820, 7830, 7840, 7860, 7865	50/case
SM4503NP	Neonatal Temp Probes for Ohmeda® IC, GC, and CP Incubators. Infant Warmer Systems 3000, 3300, and 5000	50/case
SM4504NP	Neonatal Temp Probes for GE® Giraffe OmniBed, Giraffe Incubators, Giraffe Warmers	50/case
SM4505NP	Neonatal Temp Probes for Dräger® Babytherm® 8004 & Babytherm®. 8010 CALEO Incubators	50/case
SM4506NP	Neonatal Temp Probes for Dräger® Isolette® C2000 & 8000 Incubators. Versalet 7700 Care Center	50/case
SM4507NP	Neonatal Temp Probes for ATOM Infa Warmer 1, ATOM Medical Incubators	50/case



Reflectors



SM4501NP



SM4503NP



SM4504NP



SM4505NP



SM4506NP



SM4507NP

Tympanic Temperature Probe

FEATURES AND BENEFITS

- Provides safe and secure seal of ear canal for canal to equilibrate to tympanic temperature
- Probe features a dual memory foam tip for protection of the tympanic membrane with less risk of perforation or damage

Item #	Description	Qty
SM400TY	Adult Tympanic Temperature Sensor 400 Series	50/box
SM400TYP	Pediatric Tympanic Temperature Sensor 400 Series	50/box



Temperature Probe Cable

Item #	Description	Qty
SMC400-MP-M	Standard 1/4" Phone plug 400 Series, 10ft	10/box
SMC400-P-M	Philips/HP Cable 400 Series, 10ft	10/box





SourceMark™
Medical