



SourceMark™
Medical



ANESTHESIA / AIRWAY MANAGEMENT



Enhanced Comfort Endotracheal Tube



FEATURES AND BENEFITS

- DEHP Free
- High volume, low profile cylindrical cuff
- Made of softer PVC material, reducing risk of trauma

TUBE SIZE	ITEM #	QTY
5mm	M0450CS	10/box
5.5mm	M0455CS	10/box
6mm	M0460CS	10/box
6.5mm	M0465CS	10/box
7mm	M0470CS	10/box
7.5mm	M0475CS	10/box
8mm	M0480CS	10/box
8.5mm	M0485CS	10/box
9mm	M0490CS	10/box

SuctionPlus ® Endotracheal Tube



FEATURES AND BENEFITS

- Latex & DEHP Free
- Removable standard 15mm connector
- Thermosmoothing of the tube tips
- High volume, low pressure cylindrical cuff
- One-way inflation valve with inflation pilot balloon
- Sizes 6.0-9.0mm, individual sterile package

Intended Use

SuctionPlus Endotracheal tube is designed with the evacuation lumen in order to provide protection against risk of aspiration and decrease the rate of ventilation-associated pneumonia (VAP) in the hospital.


TUBE SIZE	ITEM #	QTY
6mm	M0495C	10/box
6.5mm	M0500C	10/box
7mm	M0505C	10/box
7.5mm	M0510C	10/box
8mm	M0515C	10/box
8.5mm	M0520C	10/box
9mm	M0525C	10/box



Standard Hi-Lo Cuffed ET Tube

10 ea per box / 10 boxes per case

Hi/Lo Straight HVLP Endotracheal Tube, Cuffed, Murphy, 15mm connector



Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
ITEM #	M0430C	M0435C	M0440C	M0445C	M0450C	M0455C	M0460C	M0465C	M0470C	M0475C	M0480C	M0485C	M0490C

Oral Cuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Oral Preformed, Cuffed, Murphy, 15mm connector



Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5
ITEM #	M0530C	M0535C	M0540C	M0545C	M0550C	M0555C	M0560C	M0565C	M0570C	M0575C	M0580C	M0585C

Nasal Cuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Nasal Preformed, Cuffed, Murphy, 15mm connector



Tube Size	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0
ITEM #	M0640C	M0645C	M0650C	M0655C	M0660C	M0665C	M0670C	M0675C	M0680C

Reinforced Cuffed ET Tube

5 ea per box / 10 boxes per case

Endotracheal Tube, Reinforced, Cuffed, Murphy, 15mm connector



Tube Size	5.5	6.0	6.5	7.0	7.5	8.0	8.5
ITEM #	M0755C	M0760C	M0765C	M0770C	M0775C	M0780C	M0785C

The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

CUFFED:

- High-volume low-pressure cylindrical cuff distributes even pressure on the tracheal wall reducing chance of tissue damage
- Large cuff surface provides better seal
- Durable construction to ensure the reliability of the tube and cuff during procedures

TUBE:

- Thermo-smoothing of the tube tips to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized



Standard Uncuffed ET Tube

10 ea per box / 10 boxes per case

Straight Endotracheal Tube, Uncuffed, Murphy, 15mm connector

Tube Size	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0
ITEM #	M0420U	M0425U	M0430U	M0435U	M0440U	M0445U	M0450U	M0455U	M0460U



Oral Uncuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Oral Preformed, Uncuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0
ITEM #	M0530U	M0535U	M0540U	M0545U	M0550U	M0555U	M0560U



Nasal Uncuffed Preformed ET Tube

10 ea per box / 10 boxes per case

Endotracheal Tube, Nasal Preformed, Uncuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0
ITEM #	M0630U	M0635U	M0640U	M0645U	M0650U	M0655U	M0660U	M0665U	M0670U



Reinforced Uncuffed ET Tube

5 ea per box / 10 boxes per case

Endotracheal Tube, Reinforced, Cuffed, Murphy, 15mm connector

Tube Size	3.0	3.5	4.0	4.5	5.0
ITEM #	M0730U	M0735U	M0740U	M0745U	M0750U



The SourceMark Endotracheal Tube is a single-use medical device intended for airway management by oral and nasal intubation during procedures.

TUBE:

- Thermo-smoothing of the tube tips to reduce risk of tissue trauma
- Murphy tip and eye
- Standard tube length and depth markings
- Removable standard 15mm connector
- Fully embedded x-ray opaque strip over full length of tube

GENERAL:

- FDA 510(K) clearance
- One-way valve with inflation pilot balloon
- Sterilized




Standard Hi-Lo Cuffed ET Tube with Pre-Loaded Stylet

10 ea per box / 10 boxes per case

Tubes, Endo, Oral and Nasal, Murphy, Cuffed + Stylet, Aluminum

Tube Size	6.5	7.0	7.5	8.0
ITEM #	MX0465C	MX0470C	MX0475C	MX0480C

SourceMark's Hi-Lo cuffed Endotracheal Tube with Stylet preloaded for quick intubation



Aluminum Stylet

20 ea per box / 20 boxes per case

Indications for usage:

- Airway management during procedures
- Anesthesia

FEATURES AND BENEFITS

- Malleable aluminum which allows the stylet to be easily adapted to any desired shape
- Coated to reduce friction between the stylet and the tracheal tube to allow for easy insertion and withdrawal

Size	6FR	8FR	10FR	12FR	14FR
ITEM #	M0906	M0908	M0910	M0912	M0914



Nasopharyngeal Airways

10 ea per box / 10 boxes per case

FEATURES AND BENEFITS

- Maximum kink resistance
- Smoothly rounded edges for patient comfort
- Made of medical grade PVC
- Angled design to reduce risk of injury
- Beveled tip for easy insertion
- Available in full range of sizes

Size	14FR	16FR	18FR	20FR	22FR	24FR	26FR	28FR	30FR	32FR	34FR	36FR
ITEM #	M1014	M1016	M1018	M1020	M1022	M1024	M1026	M1028	M1030	M1032	M1034	M1036





SourceMarkTM
Medical